

## Measurement of hardness and pellet size of spherical samples

The hardness and pellet size of granulation characteristics of spherical samples such as carbon black can be measured simultaneously.

---

## Continuous measurement of different types of samples

The AS2100 can be equipped with multiple sample input ports to continuously measure samples of different equivalent sizes.

---

## Two data management

Two types of data can be saved: crushing load x pellet size distribution chart and individual load curve data.

---

Specification	
Item	Specifications
External dimensions	W355 x D560 x H340 (mm)
Weight	about 28kg
Maximum load	200cN ≒200gf
Compression speed	0.125mm/sec, 5-step variable
Pellet size	spherical 1.0≒1.7mm
Output	USB Type-B, Data transfer to PC
Power supply	AC100V±5%, 50/60Hz
Supported OS	Windows 10